



PATIENT ORDER FORM

Please complete and mail this form with your payment to:

Timm Medical Technologies, Inc.
6585 City West Parkway
Eden Prairie, MN 55344
USA

Completed by patient:

I WOULD LIKE TO ORDER THE ACTICUF™ DISPOSABLE POUCH

NAME _____

ADDRESS _____

CITY _____

STATE _____ ZIP CODE _____

TELEPHONE _____

E-MAIL _____

Check if you would like to receive Timm Medical e-mail updates

- 1 10-pack @ \$45.00 Each (\$7.50 shipping and handling)
- 2 10-packs @ \$42.00 Each (\$7.50 shipping and handling)
- 3 10-packs @ \$39.00 Each (\$7.50 shipping and handling)
- 4 or more 10-packs @ \$35.00 Each (\$7.50 shipping and handling)

Total (10-pack) Qty. = _____ x Price (see above) _____ = _____

Add Applicable Sales Tax = _____

+ \$7.50 Shipping and handling = _____

Total Amount of Order = _____

METHOD OF PAYMENT (Please do not send cash)

Check-payable to Timm Medical Technologies, Inc.

Money Order/Cashier's Check

Credit Card: Discover Visa Master Card Amex

CARD # _____ EXP. DATE _____

SIGNATURE _____ DATE _____

No refunds or exchanges for components.
Prices subject to change without notice.

QUESTIONS?
www.timmmedical.com
1-800-438-8592

ActiCuf™ Insight: Frequently Asked Questions

Q: Will people notice I'm wearing the ActiCuf product?

A: Rest assured, the ActiCuf product was engineered and carefully tested with this very concern in mind. Its unique contour forms naturally to your body. Plus, it's small enough to carry with you confidently — even in your pocket!

Q: How long can I wear the ActiCuf product?

A: One ActiCuf Disposable Pouch can be used all day however, it will need to be repositioned every three to four hours. Though rare, instances of skin irritation may occur with prolonged exposure to urine. The ActiCuf Disposable Pouch should not be worn while sleeping.

Q: Can I use a public urinal with ActiCuf product?

A: Absolutely! Because the ActiCuf product is so discreet, you can take it off and put it on again quickly and effortlessly — without anyone noticing.

Q: How can I purchase the ActiCuf product?

A: • Log on to www.timmmedical.com
• Simply call toll-free 1-800-438-8592
NOTE: Individuals with urinary incontinence should always visit their physician for a thorough evaluation.

Q: Is a prescription required?

A: No. Any male affected by incontinence who desires a renewed lifestyle is invited to experience the freedom of the ActiCuf product.

Q: Who do I contact if I have questions regarding the ActiCuf product?

A: Call toll-free 1-800-438-8592 to speak directly with an ActiCuf Professional Customer Service Representative. More product details and instructions for use can also be found on www.timmmedical.com

Distributed by:
Timm Medical Technologies, Inc.
6585 City West Parkway
Eden Prairie, MN 55344
USA

1-800-438-8592
Phone: 952-947-9410
Fax: 952-947-9411
www.timmmedical.com



PM07-0011 Rev. A 3/07
ActiCuf is a trademark of GT Urological



**Helping
you get back
to living your life**



**ACTIVE LIFESTYLE MALE
INCONTINENCE CONTROL**



The breakthrough ActiCuf™ product is a disposable, undetectable pouch designed to help men treat and manage light to moderate urinary incontinence. When placed on the penis, the padded closure gently presses down on the urethra to control urinary flow. Any leakage is immediately absorbed by the ActiCuf product, so you can take charge of your day and enjoy life on your terms...it's a new way of life.

The ActiCuf™ Disposable Pouch is so discreet that it's virtually undetectable under clothing. And its sleek, portable design makes using public urinals a worry-free experience. Removing and replacing it is fast and easy.

QUESTIONS?

www.timmmedical.com
1-800-438-8592



PROTECTIVE
ABSORBENT SHIELD

UNIQUE SHAPE
AND SIZE IS DISCREET AND
EASY TO REMOVE

PADDED CLOSURE
COMFORTABLY PREVENTS LEAKAGE

HERE'S HOW TO USE:

- 1 Hold the ActiCuf product with your thumb and forefinger.
- 2 Squeeze the ActiCuf product until the end is completely open and round in shape. Insert penis as far as possible into the opening of the device.
- 3 Release pressure so the ActiCuf product fits snugly around the penis. Adjust it as necessary to achieve the greatest comfort and ultimate leakage protection.
- 4 Using public urinals is convenient again, simply remove by:
 - Squeezing the shaft of the penis to stop urinary flow.
 - Remove the ActiCuf product to urinate
 - Replace the ActiCuf product by repeating Steps 1-3 above.

TIP: Although the ActiCuf product is engineered for the greatest ease of use, it's not uncommon for accidents to occur. Keep an extra ActiCuf product close at hand in case of sudden leaks — just until you learn to remove and replace it.

